


Technical Data

CalStudio™

Smarter Workflows. Lower Downtime.
Built for your future.



Calibration Made Simple

Meet CalStudio™, a powerful, all-in-one software solution for calibration lab management. CalStudio™ is a simple, reliable system that easily automates workflows, manages measurements from every discipline, and quickly creates procedures—no coding required. Write once, run anywhere, and simplify your lab operations.

Everything you need to write, test, record, and report – all in one place.



Three apps in one

Write and execute procedures, manage laboratory workflows, and store calibration data all in one place. Plus, cloud-based infrastructure means you can access CalStudio from anywhere.



Drive efficiency with data

Monitor KPIs like calibration time and volume to improve lab performance. Intuitive dashboards and customizable widgets allow you to track trends with the most impact.



No coding necessary

Our built-in procedure authoring software requires no coding for easy execution. Authorized users can quickly publish new tasks that sync seamlessly to technician workflows.



Easily scale over time

CalStudio™ is designed to grow with you. Whether you manage 10 or 10,000 calibrations per month, our platform can handle the complexity and volume of your operation.



Access a robust asset library

CalStudio™ comes pre-loaded with data from 800,000 common assets, so you can fill in work orders quickly.



Software by **FLUKE**

Product group:

- Electrical
- RF
- Temperature
- Humidity
- Pressure
- Flow
- Software

Streamline calibration with CalStudio

CalStudio™ is a cloud-based application for calibration lab management. Designed by the experts at Fluke Calibration, it combines three applications into one intuitive platform: a laboratory information management system (LIMS), procedure authoring and execution software, and a workflow and integration manager. Designed for end-to-end performance, CalStudio™ makes it easy for calibration labs to manage customer requests, write procedures, and complete jobs with confidence.



Customer & Asset Management

- Support for external customers/clients
- Bulk asset import from CSV file



MET/CAL Integration

- Upload and launch MET/CAL procedures
- Connect results via MET/CONNECT



Product Catalog

- Leverage Cubyt data
- Consistency & ease of use



Easy Procedure Creation

- Drag and drop
- Low-code/No-code




Author Once, Run Anywhere

- Reference agnostic
- Configuration agnostic



OptiCal™ AI Data Capture

- AI-powered reading capture
- Hands-free, no manual entry



App Integration Library

- Eliminate data segregation
- 1000's of data connectors



Custom Workflows

- Drag & drop design
- Tailored automation

Unlock Efficiency with Smart Calibration Automation



2 Days
Reduced TAT

Faster turnaround,
quicker results.



+\$14,000
Invoiced / Tech

Drive higher sales
per tech / year.



+2 to 5%
Margin Increase

Improve your
overall bottom line.



800k Devices
Preloaded

Minimal
onboarding time.

Fluke. Keeping your world up and running.™

[Fluke.com](https://www.fluke.com)

©2026 Fluke Corporation.
Specifications subject to change without notice.
260373-en

Modification of this document is not permitted
without written permission from Fluke Corporation.

FLUKE